



11

Aspiration, Sonden

Aspiration, Probes

Aspiration, Sondes

Aspiración, Sondas

Aspirazione, Sonde



Aspiration, Sonden Aspiration, Probes



Trokare
Trocars
Trocarts
Trócares
Trequarti



1:1

1:1

UNIVERSAL

	Fig. 1	Fig. 2	Fig. 3	Fig. 4
	6.0 mm Ø	4.5 mm Ø	3.0 mm Ø	1.75 mm Ø
38-230-003		x	x	x
38-230-004	x	x	x	x

Satz
Set
Jeu
Juego
Serie



1:1

1:1

1:1

1:1

UCKERMANN

	8 mm Ø	10 mm Ø
130 mm – 5"	38-380-008	38-380-010

1:2



SINUS

	3.5 mm Ø	5.0 mm Ø	5.0 mm Ø
170 mm – 6 3/4"	38-384-035	38-386-050	38-390-033

1:2



1:1

1:1

1:1



Aspiration, Sonden Aspiration, Probes

Trokare
Trocars
Trocarts
Trócares
Trequarti




LICHTWITZ

	1.5 mm Ø	1.8 mm Ø	2.0 mm Ø
190 mm – 7 1/2"	346-336-150	346-366-180	346-366-200

1:2

OCHSNER


		100 mm – 4"
4 mm Ø	12 Char./Fr.	38-282-012
4 2/3 mm Ø	14 Char./Fr.	38-282-014
5 1/3 mm Ø	16 Char./Fr.	38-282-016
6 mm Ø	18 Char./Fr.	38-282-018
6 2/3 mm Ø	20 Char./Fr.	38-282-020
7 1/3 mm Ø	22 Char./Fr.	38-282-022
8 mm Ø	24 Char./Fr.	38-282-024
8 2/3 mm Ø	26 Char./Fr.	38-282-026



NELSON

	7.5 mm Ø	10 mm Ø	11 mm Ø	11.6 mm Ø
210 mm – 8 1/4"	322-166-250	322-166-300	322-166-330	322-166-350


1:2



Aspiration, Sonden Aspiration, Probes



Trokare
Trocars
Trocart
Trócares
Trequarti



LANDAU

1:2

	2 mm Ø
160 mm – 6 ¼"	418-466-160



KRAUSE

1:2

	5.5 mm Ø
150 mm – 6"	418-501-150



PIERCE

	170 mm – 6 ¾"
3 mm Ø	322-119-300
4 mm Ø	322-119-400
5 mm Ø	322-119-500



COAKLEY

1:2

	2.25 mm Ø
160 mm – 6 ¼"	322-136-020

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Saugrohre
Suction tubes
Canules d'aspiration
Tubos de aspiración
Tubi per aspirazione



COOLEY

	8 mm Ø		10 mm Ø		6 mm Ø		8 mm Ø
350 mm – 13 ¾"	418-557-350	310 mm – 12 ¼"	418-559-310	300 mm – 12"	418-561-300	330 mm – 13"	418-563-330

Aspiration, Sonden Aspiration, Probes



Saugrohre
Suction tubes

Canules d'aspiration
Tubos de aspiración
Tubi per aspirazione



POOLE

1:3

		220 mm – 8"
8 mm Ø	23 Charr./Fr.	38-551-008



POOLE

1:3

		220 mm – 8"
10 mm Ø	30 Charr./Fr.	38-550-010



ANDREWS-PYNCHON

1:3

235 mm – 9 1/4"	38-523-235
-----------------	------------



DE BAKEY

1:3

275 mm – 10 3/4"	304-150-270
------------------	-------------



38-521-270/2



38-521-270/1



YANKAUER

1:3

270 mm – 10 3/4"	38-521-270
------------------	------------

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Ohrkatheter

Eustachian catheters
Sondes auriculaires
Sondas auriculares
Sonde auricolari



HARTMANN

1:2

1:2

	Fig. 1	Fig. 2
2.0 mm Ø	437-797-200	437-798-200
2.5 mm Ø	437-797-250	437-798-250
3.0 mm Ø	437-797-300	437-798-300



EICKEN-KILLIAN

1:2

130 mm – 6"	2.0 mm Ø	38-685-020
130 mm – 6"	2.5 mm Ø	38-685-025
130 mm – 6"	3.0 mm Ø	38-685-030
130 mm – 6"	4.0 mm Ø	38-685-040



HARTMANN

1:1

1:1

1:1

1:1

1:1


1:1

	Fig. 0	Fig. 1	Fig. 2	Fig. 3	Fig. 4	Fig. 5
150 mm – 6"	38-621-000	38-621-001	38-621-002	38-621-003	38-621-004	38-621-005

Aspiration, Sonden Aspiration, Probes



Ohrkatheter
Eustachian catheters
Sondes auriculaires
Sondas auriculares
Sonde auricolari



436-336-000

Schlaucholive für Ohrkatheter
Tube olive for eustachian catheter
Olive pour sonde auriculaire
Olive para sondas auriculares
Olive per sonde auricolari

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Ohrspritze
Ear syringe
Seringe auriculaire
Jeringa auricular
Siringhe auricolari



50 ccm 2 oz	44-310-050
75 ccm 3 oz	44-310-075
100 ccm 4 oz	44-310-100
150 ccm 6 oz	44-310-150
200 ccm 8 oz	44-310-200

komplett mit 2 Aufsätzen
complete with 2 nozzles
complet avec 2 canules
completa con 2 cánulas
completa con 2 terminali



44-310-002



44-310-005




44-310-008

Aspiration, Sonden Aspiration, Probes



Saug- und Spülkanülen
Suction and irrigation cannulae
Canules pour aspiration et lavage
Cánulas para aspiración y lavado
Cannule per aspirazione e lavaggio



EICKEN

2.5 mm Ø	38-687-025	38-689-025
3.0 mm Ø	38-687-030	38-689-030
4.0 mm Ø	38-687-040	38-689-040

leicht gebogen
slightly curved
légèrement courbé
ligeramente curvada
leggermente curva

stark gebogen
strongly curved
très courbé
muy curvada
fortemente curva



LUER-LOCK

3.0 mm Ø
38-593-030

LUER-LOCK

2.5 mm Ø
38-591-030

FISCH

1.8 mm Ø	38-590-018
2.5 mm Ø	38-590-025
3.0 mm Ø	38-590-030
4.5 mm Ø	38-590-040

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Saug- und Spülkanülen
Suction and irrigation cannulae
Canules pour aspiration et lavage
Cánulas para aspiración y lavado
Cannule per aspirazione e lavaggio



0.8 mm Ø	350-288-000	350-289-900	0.8 mm Ø	350-289-000	350-289-450
----------	-------------	-------------	----------	-------------	-------------

Tonsillenkanüle
Tonsil needle
Aiguille pour amygdales
Aguja para amígdalas
Cannule per tonsille

gerade
straight
droite
recta

gebogen
curved
courbe
curva

Septumkanüle
Septum needle
Aiguille pour le cloison nasal
Aguja para el tabique nasal
Ago per setto nasale

bajonettförmig
bayonet
baionette
bayonetta

abgewinkelt
angled on flat
coudé sur le plat
angulado
angolato

Aspiration, Sonden Aspiration, Probes



11

Saugrohre
Suction tubes
Tubes d'aspiration
Tubos para aspiración
Cannule d'aspirazione



ZÖLLNER

0.7 mm Ø
344-095-000

komplett mit 1 aufsetzbaren Kanüle, 0.7 mm Ø
complete with 1 slip-on cannula, 0.7 mm Ø
complet avec 1 canule démontable, 0.7 mm de Ø
completo con 1 cánula desmontable de 0.7 mm Ø
completa di 1 cannula applicabile di mm 0.7 di Ø




MAGILL

	245 mm – 9 3/4"
3.0 mm Ø	304-151-300
3.5 mm Ø	304-151-350
4.0 mm Ø	304-151-400

IOWA

1:2

	0.7 mm Ø
160 mm – 6 1/4"	344-090-160



ROSEN

1:2

1:2

	x = 80 mm	x = 65 mm
0.5 mm Ø	38-572-005	38-573-005
0.6 mm Ø	38-572-006	38-573-006
0.7 mm Ø	38-572-007	38-573-007
0.9 mm Ø	38-572-009	38-573-009
1.0 mm Ø	38-572-010	38-573-010
1.2 mm Ø	38-572-012	38-573-012
1.4 mm Ø	38-572-014	38-573-014
1.5 mm Ø	38-572-015	38-573-015
1.8 mm Ø	38-572-018	38-573-018
2.0 mm Ø	38-572-020	38-573-020
2.5 mm Ø	38-572-025	38-573-025
3.0 mm Ø	38-572-030	38-573-030

allgaier
instrumente
medical


Aspiration, Sonden Aspiration, Probes

Saugrohre
Suction tubes
Tubes d'aspiration
Tubos para aspiración
Cannule d'aspirazione



BARON

		75 mm – 3"	100 mm – 4"
1 mm Ø	3 Charr./Fr.	344-097-030	344-097-030/L
1 1/2 mm Ø	5 Charr./Fr.	344-097-050	344-097-050/L
2 1/2 mm Ø	7 Charr./Fr.	344-097-070	344-097-070/L



HOUSE

	90 mm – 3 1/2"
0.9 x 0.7 mm Ø	344-099-010
1.2 x 0.9 mm Ø	344-099-020
1.5 x 1.2 mm Ø	344-099-030
1.8 x 1.8 mm Ø	344-099-040
2.5 x 2.0 mm Ø	344-099-050
3.0 x 2.0 mm Ø	344-099-060
4.5 x 2.0 mm Ø	344-099-070



FERGUSSON

		170 mm – 6 3/4"
2 mm Ø	6 Charr./Fr.	418-515-020
2 1/2 mm Ø	8 Charr./Fr.	418-515-025
3 mm Ø	9 Charr./Fr.	418-515-030
3 1/2 mm Ø	10 Charr./Fr.	418-515-035
4 mm Ø	12 Charr./Fr.	418-515-040
5 mm Ø	15 Charr./Fr.	418-515-050


1:2

Aspiration, Sonden Aspiration, Probes



11

Saugrohre
Suction tubes
Tubes d'aspiration
Tubos para aspiración
Cannule d'aspiratione



PLESTER

		170 mm – 6 3/4"
1.6 mm Ø	5 Charr./Fr.	418-517-015
2 mm Ø	6 Charr./Fr.	418-517-020
2 1/2 mm Ø	8 Charr./Fr.	418-517-025
3 mm Ø	9 Charr./Fr.	418-517-030

FRAZIER

		170 mm – 6 3/4"
2 mm Ø	6 Charr./Fr.	38-571-006
2 1/2 mm Ø	7 Charr./Fr.	340-291-070
2 3/4 mm Ø	8 Charr./Fr.	38-571-008
3 mm Ø	9 Charr./Fr.	340-291-090
3 1/2 mm Ø	10 Charr./Fr.	38-571-010
4 mm Ø	12 Charr./Fr.	38-571-012
5 mm Ø	15 Charr./Fr.	340-291-150



ADSON

		140 mm – 5 1/4"	210 mm – 8 1/4"
3 mm Ø	9 Charr./Fr.	340-286-140	
4 mm Ø	12 Charr./Fr.		340-286-210
5 mm Ø	15 Charr./Fr.		340-286-211

200 mm – 8"	437-835-200
-------------	-------------

Saug-Raspatorium – biegsam, scharf
Suction raspatorium – malleable, sharp
Rugine aspirante – flexible, tranchant
Legra aspiradora – maleable y agudo
Elevatore-aspiratore – flessibile, tagliente

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Redon Führungsnadeln
Guide needles
Aiguilles d'introduction
Agujas de guía
Aghi guida




		195 mm – 5 3/4"
2 1/2 mm Ø	8 Charr./Fr.	18-581-008
3 1/2 mm Ø	10 Charr./Fr.	18-581-010
4 mm Ø	12 Charr./Fr.	18-581-012
4 1/2 mm Ø	14 Charr./Fr.	18-581-014
5 1/2 mm Ø	16 Charr./Fr.	18-581-016
6 mm Ø	18 Charr./Fr.	18-581-018



		195 mm – 5 3/4"
2 1/2 mm Ø	8 Charr./Fr.	18-585-008
3 1/2 mm Ø	10 Charr./Fr.	18-585-010
4 mm Ø	12 Charr./Fr.	18-585-012
4 1/2 mm Ø	14 Charr./Fr.	18-585-014
5 1/2 mm Ø	16 Charr./Fr.	18-585-016
6 mm Ø	18 Charr./Fr.	18-585-018

Lanzettenspitze
Lancet tip
Pointe lancéolaire
Punta lanceolada
Punta lanceolata



		195 mm – 5 3/4"
2 1/2 mm Ø	8 Charr./Fr.	18-583-008
3 1/2 mm Ø	10 Charr./Fr.	18-583-010
4 mm Ø	12 Charr./Fr.	18-583-012
4 1/2 mm Ø	14 Charr./Fr.	18-583-014
5 1/2 mm Ø	16 Charr./Fr.	18-583-016
6 mm Ø	18 Charr./Fr.	18-583-018

Aspiration, Sonden Aspiration, Probes



Injektionskanülen
Hypodermic needles

Aiguilles hypodermiques
Agujas hipodérmicas
Aghi hypodermici



PRAVAZ (LUER-LOCK)

Fig. 1	0.90 x 38 mm Ø	418-598-010
Fig. 2	0.80 x 32 mm Ø	418-598-020
Fig. 12	0.70 x 32 mm Ø	418-598-120
Fig. 14	0.65 x 28 mm Ø	418-598-140
Fig. 16	0.60 x 25 mm Ø	418-598-160
Fig. 18	0.50 x 23 mm Ø	418-598-180
Fig. 20	0.45 x 20 mm Ø	418-598-200



Packung à 12 Stück
Package of 12 pieces
Paquet de 12 pièces
Paquete de 12 piezas
Confezione da 12 pezzi

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Biopsiekanülen
Biopsy needles

Canules pour biopsie
Agujas para biopsia
Aghi per biopsie

QUINKE

0.70 x 50.0 mm Ø	322-230-070
0.70 x 76.0 mm Ø	322-226-070
0.70 x 89.0 mm Ø	322-228-070
1.00 x 76.0 mm Ø	322-226-100
1.00 x 89.0 mm Ø	322-228-100
1.20 x 76.0 mm Ø	322-226-120
1.20 x 89.0 mm Ø	322-228-120

TUOHY

1.20 x 76.0 mm Ø	418-628-120
1.40 x 76.0 mm Ø	418-628-140
1.60 x 76.0 mm Ø	418-628-160

BIER

1.00 x 80.0 mm Ø	322-232-100
1.20 x 80.0 mm Ø	322-232-120
1.50 x 80.0 mm Ø	322-232-150
1.80 x 80.0 mm Ø	322-232-180

KLIMA-ROSEGGER

1.80 x 35.0 mm Ø	418-656-180
------------------	-------------

VERRES

2.00 x 100 mm Ø	92-010-100
2.00 x 120 mm Ø	92-010-120
2.00 x 130 mm Ø	92-010-130
2.00 x 140 mm Ø	92-010-140
2.00 x 150 mm Ø	92-010-150

Aspiration, Sonden Aspiration, Probes



11

Biopsiekanülen Biopsy needles

Canules pour biopsie
Agujas para biopsia
Aghi per biopsie



MENGHINI

1.00 x 35.0 mm Ø	38-420-001
1.00 x 70.0 mm Ø	38-420-002
1.20 x 70.0 mm Ø	38-420-003
1.20 x 100.0 mm Ø	38-420-004
1.40 x 70.0 mm Ø	38-420-005
1.40 x 100.0 mm Ø	38-420-006
1.60 x 70.0 mm Ø	38-420-007
1.80 x 100.0 mm Ø	38-420-008



MENGHINI

1.40 x 168.0 mm Ø	38-430-140
1.80 x 168.0 mm Ø	38-430-180



Aspiration, Sonden Aspiration, Probes

Verbindungsstücke und Hähne Connecting pieces and stopcocks

Raccords et robinets
Piezas de conexión y llaves
Raccordi e rubinetti



1:1

418-836-000



1:1

418-841-000

Aspiration, Sonden Aspiration, Probes



Biopsiekanülen Biopsy needles

Canules pour biopsie
Agujas para biopsia
Aghi per biopsie



Außendurchmesser
Outer diameter
Diamètre externe
Diámetro externo
Diametro esterno



Innendurchmesser
Inner diameter
Diamètre interne
Diámetro interno
Diametro interno



FRANKLIN-SILVERMAN

⊗	⊙	
2.00 x 105.0 mm Ø	1.50 x 155.0 mm Ø	322-326-200
2.00 x 130.0 mm Ø	1.50 x 180.0 mm Ø	322-326-140

Zur perkutanen Nierenbiopsie, einschl. Suchnadel
For kidney biopsy, incl. guide needle
Pour la biopsie des reins, aiguille d'introduction incl.
Para la biopsia de los riñones, con aguja de guía
Per la biopsia renale, con ago guida

Aspiration, Sonden Aspiration, Probes

Biopsiekanülen Biopsy needles

Canules pour biopsie
Agujas para biopsia
Aghi per biopsie



Außendurchmesser
Outer diameter
Diamètre externe
Diámetro externo
Diametro esterno



Innendurchmesser
Inner diameter
Diamètre interne
Diámetro interno
Diametro interno



FRANKLIN-SILVERMAN

⊗	⊙		
2.00 x 100.0 mm Ø	1.50 x 150.0 mm Ø	38-450-020	38-454-020
2.50 x 145.0 mm Ø	1.50 x 195.0 mm Ø	38-450-025	38-454-025

Zur Leber- und Nierenbiopsie
For liver and kidney biopsy
Pour la biopsie des reins et de la foie
Para la biopsia de los riñones y del hígado
Per la biopsia renale ed epatica

ohne Widerhaken
without barb
sans crochet
sin garfio
senza gancio

mit Widerhaken
barbed
avec crochet
con garfio
con gancio

⊗	⊙		
2.00 x 100.0 mm Ø	1.50 x 150.0 mm Ø	322-348-150	322-352-150
2.50 x 145.0 mm Ø	1.50 x 195.0 mm Ø	322-348-200	322-352-200

Kanüle, allein
Biopsy needle, only
Aiguille pour biopsie, seul
Aguja para biopsia, sola
Ago per biopsia, solo

Aspiration, Sonden Aspiration, Probes



Glas/Metall Spritzen
Glas/metal syringes

Seringues en verre/métal
Jeringas de vidrio/metal
Siringhe vetro/metallo



Ersatzteile und Zubehör auf Anfrage
Spare parts and accessories upon request
Pièces de rechange et accessoires sur demande
Piezas de cambio y accesorios según demanda
Pezzi di ricambio ed accessori a richiesta



LUER

1 ml	418-896-010	418-896-010/1
2 ml	418-896-020	418-896-020/1
5 ml	418-896-050	418-896-050/1
10 ml	418-896-100	418-896-100/1
20 ml	418-896-200	418-896-200/1
50 ml	418-896-500	418-896-500/1

Spritze, komplett
Syringe, complete
Seringue, complète
Jeringa, completa
Siringa, completa

Zylinder
Barrel
Cylindre
Cilindro
Cilindro



LUER-LOCK

1 ml	418-906-010	418-906-010/1
2 ml	418-906-020	418-906-020/1
5 ml	418-906-050	418-906-050/1
10 ml	418-906-100	418-906-100/1
20 ml	418-906-200	418-906-200/1
50 ml	418-906-500	418-906-500/1

Spritze, komplett
Syringe, complete
Seringue, complète
Jeringa, completa
Siringa, completa

Zylinder
Barrel
Cylindre
Cilindro
Cilindro



LUER-LOCK

2 ml	350-286-020	350-286-021
5 ml	350-286-050	350-286-051
10 ml	350-286-100	350-286-101
20 ml	350-286-200	350-286-201

Spritze, komplett
Syringe, complete
Seringue, complète
Jeringa, completa
Siringa, completa

Zylinder
Barrel
Cylindre
Cilindro
Cilindro

allgaier
instrumente
medical

Aspiration, Sonden Aspiration, Probes

Wund- und Blasenspritzen
Wound and bladder syringes

Seringues pour plaies et vessie
Jeringas para enjuagar vejiga y heridas
Siringhe per ferite e vescica



1:2
418-926-004



418-926-002

JANET

50 ml	418-926-050	418-926-050/1
75 ml	418-926-075	418-926-075/1
100 ml	418-926-100	418-926-100/1
150 ml	418-926-150	418-926-150/1
200 ml	418-926-200	418-926-200/1
250 ml	418-926-250	418-926-250/1

Spritze, komplett
Syringe, complete
Seringue, complète
Jeringa, completa
Siringa, completa

Zylinder, allein
Barrel, only
Cylindre, seul
Cilindro, sola
Cilindro, solo

Aspiration, Sonden Aspiration, Probes



Sonden
Probes and directors
Sondes et stylets
Sondas y estiletes
Specilli e sonde



	1.0 mm Ø	1.5 mm Ø	2.0 mm Ø
115 mm – 4 ½"	33-118-115	33-119-115	33-120-115
130 mm – 5"	33-118-130	33-119-130	33-120-130
145 mm – 5 ¾"	33-118-145	33-119-145	33-120-145
160 mm – 6 ¼"	33-118-160	33-119-160	33-120-145
180 mm – 7"	33-118-180	33-119-180	33-120-180
200 mm – 8"	33-118-200	33-119-200	33-120-200
250 mm – 9 ¾"	33-118-250	33-119-250	33-120-250
300 mm – 12"	33-118-300	33-119-300	33-120-300



2 mm Ø	
115 mm – 4 ½"	33-140-115
130 mm – 5"	33-140-130
145 mm – 5 ¾"	33-140-145
160 mm – 6 ¼"	33-140-160
180 mm – 7"	33-140-180
200 mm – 8"	33-140-200
250 mm – 9 ¾"	33-140-250



2 mm Ø	
115 mm – 4 ½"	33-170-120
130 mm – 5"	33-170-130
145 mm – 5 ¾"	33-170-145
160 mm – 6 ¼"	33-170-160

Myrtenblattsonde
myrtle-shaped
en forme de myrte
en forma de mirto
forma di mirto

mit Öhr
with eyelet
avec oeil
con ojo
con cruna



115 mm – 4 ½"	33-180-115
130 mm – 5"	33-180-130
145 mm – 5 ¾"	33-180-145
160 mm – 6 ¼"	33-180-160
180 mm – 7"	33-180-180
200 mm – 8"	33-180-200



115 mm – 4 ½"	33-182-115
130 mm – 5"	33-182-130
145 mm – 5 ¾"	33-182-145
160 mm – 6 ¼"	33-182-160
180 mm – 7"	33-182-180
200 mm – 8"	33-182-200



NÉLATON 1:2

160 mm – 6 ½"	
33-183-160	




DOYEN

130 mm – 5"	33-187-130
145 mm – 5 ¾"	33-187-145
160 mm – 6 ¼"	33-187-160

allgaier
instrumente
medical


Aspiration, Sonden Aspiration, Probes

Watteträger
Cotton applicators
Porte-cotons
Porta-algodones
Portacotoni




FARRELL-LATHBURY

120 mm – 4 ¾"	0.9 mm Ø	33-520-120	320-288-120	33-530-120
140 mm – 5 ½"	0.9 mm Ø	33-520-140	320-288-140	33-530-140
165 mm – 6 ¼"	1.1 mm Ø	33-520-165	320-288-160	33-530-165
180 mm – 7"	1.2 mm Ø	33-520-180	320-288-180	33-530-180
230 mm – 9"	1.3 mm Ø	33-520-230	320-288-230	33-530-230



	1.2 mm Ø	1.5 mm Ø	1.5 mm Ø
120 mm – 4 ¾"	416-406-120	320-286-120	320-288-120
140 mm – 5 ½"	416-406-140	320-286-140	320-288-140
160 mm – 6 ¼"	416-406-160	320-286-160	320-288-160
180 mm – 7"	416-406-180	320-286-180	320-288-180
230 mm – 9"			320-288-230



	2.5 mm Ø	2.5 mm Ø
280 mm – 11"	416-536-280	416-538-280

Aspiration, Sonden Aspiration, Probes



Watteträger
Cotton applicators
Porte-cotons
Porta-algodones
Portacotoni



**JOBSON-
HORNE**

180 mm – 7"
33-620-180



FEIN 1:2

280 mm – 11"
33-551-280



PELTESON 1:2

230 mm – 9"
416-487-230